

Contact Us

Capital Biosciences, Inc.

900 Clopper Rd, Suite120

Gaithersburg, MD 20878

Phone: **1-800-475-2812**

Email: info@capitalbiosciences.com

Web: www.capitalbiosciences.com

SubX™ RNA Isolation Kit

INSTRUCTION MANUAL (v1.0)

Catalog No. VIR-0100

Up to 100 isolations.

For Research Use Only

Product Description

SubX™ kit is designed for isolation of circulating cell-free RNA (cfRNA) and viral RNA directly from the medium. Our technology is based on the use of proprietary bi-functional substance (SubX™) that binds DNA and RNA under physiological conditions (e.g. directly in biological liquids), followed by adsorption on a solid phase matrix. Since SubX™ captures DNA and RNA via phosphate residues (groups) it allows for elimination of bias related to both AT/CG content and fragments length, thus improving extraction efficacy and accuracy of downstream applications.

Kit components, shipping and storage conditions:

- **SubX™ Solution** (2 x 7 ml). Ready-to-use. **Store refrigerated +4°C.**
- **Binding Matrix** (2 x1 ml). Ready-to-use. Store at room temperature in a dark place.
- **Wash-1 Solution** (60 ml). Ready-to-use. Store at room temperature in a dark place.
- **Wash-2 Solution** (2 x 11 ml). **Add 55 ml of 96% molecular grade ethanol before use.** Store at room temperature in a dark place.
- **Elution Buffer:** (5 ml); Ready-to-use. **Store refrigerated +4°C.**

Kits are shipped at ambient temperatures.

Important Notes

- Optimal results are achieved when Wash-2 solution is diluted with molecular grade 96% ethanol.
- We recommend using Posi-Click 1.5 ml Eppendorf tubes (Denville Scientific, Inc., Cat. C2170) to prevent material loss during vortexing due to possible cap opening.
- Binding Matrix can be separated from the liquid phase with magnetic stand such as DynaMag™-2 Magnet. **Very important: always pellet Binding Matrix by centrifugation after step 3 of the protocol.** If using magnetic separation protocol, centrifugation for a short time is necessary to collect all materials from the tube cap.

Protocol:**RNA Isolation from 1 ml or 2 ml medium (cell culture, VTM, etc).**

(100 isolations)

- 1) Add **100 µl** of **SubX™**, to 1 ml VTM, close the tube and immediately vortex for 10-15 seconds. Turbidity may occur. Incubate 10-15 min at room temperature; vortex 3-4 times during incubation.
- 2) Vortex **Binding Matrix** for 15 sec; open the tube and add **20 µl** of **Binding Matrix** slurry to **SubX™**-medium mixture, close the tube and incubate 5-10 min at room temperature; briefly (~10 sec) vortex 3-4 times during incubation.
- 3) Centrifuge the tube for 5 min at 14,000 rpm, aspirate the supernatant carefully by pipetting taking care not to disturb the **Binding Matrix**, discard supernatant.
- 4) Add **0.4 ml Wash-1** solution to the pellet, close the tube and briefly (~10 sec) vortex and incubate for 5 min;
- 5) Add **0.2 ml 96% Ethanol**, close the tube and briefly (~10 sec) vortex and incubate for 5 min.
- 6) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 7) Add **0.6 ml Wash-2** solution to the pellet, close the tube and briefly (~10 sec) vortex and incubate for 5 min
- 8) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 9) Repeat steps 7 and 8 (once or twice)**
- 10) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 11) **ESSENTIAL step!** Briefly centrifuge the tube for 15 sec at 14,000 rpm to remove residual traces of washing buffer from the tube wells, then aspirate traces of Wash-2 solution by using a 10-20 µl micropipette tip and discard.
- 12) Add 30-50 µl of **Elution Buffer** to the pellet, tightly close the tube and incubate at +50°C for 5 min. Briefly (~10 sec) vortex the tube 5-6 times during incubation.
- 13) Centrifuge the tube for 1 min at 14,000 rpm, carefully aspirate the RNA containing supernatant, transfer to a new tube and **save**.
- 14) **(Optional)** Repeated elution can yield additional 10-20% of RNA.
- 15) Keep the RNA containing supernatant at -20°C or -80°C for long-term storage.

Protocol:**RNA Isolation from saliva.**

- 1) Centrifuge crude saliva to pellet solid food remains
- 2) Aspirate carefully supernatant not to disturb the pellet and transfer to new tube.
- 3) Add equal volume of PBS (e.g. 0.5 ml saliva + 0.5 ml PBS)
- 4) Add **140 µl** of **SubX™** solution, close the tube and briefly (~10 sec) and immediately vortex. Turbidity may occur, vortex mixture and incubate 5 min at room temperature; vortex 3-4 times during incubation.
- 5) Optional: sample can be stored at this point for at least 72 hours at room temperature before processing further, if necessary (lack of space or time).
- 6) Vortex **Binding Matrix** for 15 sec; open the tube and add **20 µl** of **Binding Matrix** slurry to **SubX™**-medium mixture, close the tube and incubate 5-10 min at room temperature; briefly (~10 sec) vortex 3-4 times during incubation.
- 7) Centrifuge the tube for 5 min at 14,000 rpm, aspirate the supernatant carefully by pipetting taking care not to disturb the **Binding Matrix**, discard supernatant.
- 8) Add **0.4 ml Wash-1** solution to the pellet, close the tube and briefly (~10 sec) vortex and incubate for 5 min;
- 9) Add **0.2 ml 96% Ethanol**, close the tube and briefly (~10 sec) vortex and incubate for 5 min.
- 10) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 11) Add **0.6 ml Wash-2** solution to the pellet, close the tube and briefly (~10 sec) vortex and incubate for 5 min
- 12) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 13) Repeat steps 7 and 8 (once or twice)**
- 14) Centrifuge the tube for 1 min at 14,000 rpm, aspirate and discard supernatant;
- 15) **ESSENTIAL step!** Briefly centrifuge the tube for 15 sec at 14,000 rpm to remove residual traces of washing buffer from the tube wells, then aspirate traces of Wash-2 solution by using a 10-20 µl micropipette tip and discard.
- 16) Add 30-50 µl of **Elution Buffer** to the pellet, tightly close the tube and incubate at +50°C for 5 min. Briefly (~10 sec) vortex the tube 5-6 times during incubation.
- 17) Centrifuge the tube for 1 min at 14,000 rpm, carefully aspirate the RNA containing supernatant, transfer to a new tube and **save**.
- 18) **(Optional)** Repeated elution can yield additional 10-20% of RNA.
- 19) Keep the RNA containing supernatant at -20°C or -80°C for long-term storage.