

Osmium Plasma Coater

For coating of conductive thin films
on scanning electron microscopy
samples



Overview




No charging. No grains. No heat damages.
An ultimate solution to obtain clear SEM images.

Osmium Plasma Coater

The Osmium Plasma Coater is a kind of plasma CVD coater which is utilizing a method for coating samples in the negative glow phase domain of DC plasma discharging. As conventional types of metal coaters, there are various kinds of so called sputter coater which use the heavy metals as coating material, such as gold, platinum, palladium etc. With these methods, the problem of granularity of coating itself is inevitable by any means. Also, if examined at high magnification, the specimens coated by these methods can not avoid charging, heat damage and contamination caused by exposure to strong electron beam with SEM, resulting in lower resolution of SEM image.

A gas reactor in which the anode and cathode plates are placed vertically is evacuated to a high vacuum, then sublimated osmium tetroxide (OsO_4) is introduced, and a DC glow discharge is generated under certain conditions, while maintaining a predetermined gas pressure. The osmium metal molecules excited by the collision of electrons instantaneously become plasma between the two electrodes. In particular, in the negative glow phase domain, there is strong light emission due to the fierce diffusion of the concentrated positively ionized metal molecules. As a result, on the surface of the specimen placed in this negative glow phase, the positive ion metal molecules uniformly adhere to the specimen surface and form a perfectly amorphous metal coating of molecular level, which molds the microstructure of the specimen surface faithfully.



Applications

Osmium thin film

- Conductive thin film for SEM specimen
- Prevention of contamination for SEM/TEM specimen
- Protective film for AFM specimen
- Conductive protective film for SPM specimen
- Protective film for SPM cantilever

Plasma-polymerized film (naphthalene)

- Protective film for FIB specimen
- Prevention for peeling sample from embedding resin
- Coating for fluorine resin (Os-PF hybrid coating)
- Support film for TEM grid
- Drift prevention for ultra-thin sections for TEM

Osmium ultra-thin film [optional]

- Observation of the ultrafine structure of insulators by FE-SEM
- Quantitative analysis of top surface of insulators by ESCA/AES
- Enhancement of conductivity for TEM specimen
- Static prevention treatment for AFM specimen
- Antistatic treatment for STM specimen
- Etching (Mixed gas method only)



Hydrophilizing treatment [optional]

- Pre-hydrophilization treatment for osmium coating on fluoro resin surface etc.
- Prevention for peeling sample from embedding resin
- Hydrophilization of the support film surface for TEM
- Hydrophilization of grids for TEM
- Hydrophilization of a diamond knife for ultramicrotome
- Improvement of wettability

Deep well electrode [optional]

- Coating for samples with a height

Multiple safety measures provide operator protection from OsO₄ exposure.



Benefits

Osmium film

- No grains:** amorphous metal coating is formed
- No heat damage:** coating is proceed under room temp.
- No electron beam damage:** the melting temperature of Osmium is 2700 deg C
- No contamination:** coating is started from vacuum state

Plasma-polymerized film (naphthalene)

- Strong film:** withstand the gallium ions beams used in FIB
- Heat resistant and insulated film**
- No grains:** amorphous metal coating is formed
- No heat damage:** coating is proceed under room temp.
- No electron beam damage**
- No contamination:** coating is started from vacuum state

Usability

- Automatic:** whole process is completed just by pressing START button after settings
- Intuitive:** coating thickness is controlled by thickness setting, not by discharge time
- Short coating time:** a few nm/a few seconds
- Easy to exchange Os ampoule:** detachable reservoir, observation window (for checking of remaining amount of Os) and built-in ampoule cutter on reservoir

Safety


- Fully automated:** reduce chance of human error
- Safety features on Os reservoir:** gas port integrated locking pin and built-ion ampoule cutter
- Interlocking system with reaction chamber:** unable to open the chamber unless OsO₄ is exhausted, unable to introduce OsO₄ when the chamber open
- Failsafe to prevent Os leak at power down**
- Osmium absorption filter:** not require any ventilator

Safety features on osmium reservoir


Filgen's OsO₄ reservoir is equipped with some safety measures including a mechanical interlocking system for safety operation.

- Gas-tight and robust design
- Safety interlocking system with gas port integrated locking pin (see below)
- Built-in ampoule cutter: can reduce risks for exposure to toxic OsO₄ gas.
- Observation window: can check remaining amount of OsO₄ crystals
- Detachable from main unit and can be stored in freezer


The gas port integrated locking pin of reservoir can simultaneously control status for reservoir locking and gas supply in safe.



Full-opened (during operation)
Reservoir: locked (und detachable) by manually inserted locking pin of reservoir
OsO₄ gas: can be introduced (depends on solenoid status)



Half-closed (before reservoir detach)
Reservoir: locked (und detachable) by manually inserted locking pin of reservoir
OsO₄ gas: blocked *residual gas in the inlet pipe must be evacuated at this step



Full-closed (after operation)
Reservoir: unlocked (detachable)
OsO₄ gas: blocked

Product line-up

Optional units

Osmium film dedicated model #OPC60A



Coating type: Osmium
Thickness: several nanometer to hundreds nanometer

Conductive ultra-thin film coating system



Coating type: Osmium
Thickness: 0.5 to 3.0 nanometer

Low-current method

Forming a film with very low discharge current by gradually introducing OsO_4 gas to highly vacuumed chamber with applying voltage. This enables to enhance reproducibility of coating thickness by lowering uncontrollable over current generated at initial phase of plasma discharge

Mixed-gas method

Forming a film slowly with relatively low concentration of OsO_4 gas by mixing an inactive gas to chamber. By discharging only with air or other inactive gas, it is also possible to etch your sample.

Osmium / Plasma-polymerized film model #OPC80T



Coating type: Osmium / PF (naphthalene)
Thickness: several nanometer to hundreds nanometer

Hydrophilizing treatment system



Enhance adhesiveness of coating to your samples by applying this feature before coating

Deep well electrode



Enable to coat your samples that have a height up to 44 mm.

Accessories

Osmium reservoir



Detachable, equipped with safety features.

Transportation Container





Enable safer transportation for OsO_4 reservoir. Stainless steel made.

OsO_4 ampoules



With exclusive optional units, Filgen's OPC provide a variety of pretreatment methods for your EM samples



Specifications – base units

Type	Osmium / Plasma-polymerized film model	Osmium film dedicated model
Model #	OPC80T	OPC60A
Reaction chamber	Glass chamber, 160(φ) x 105(H) mm	Glass chamber, 120(φ) x 73(H) mm
Max. sample size	44(W) x 44(D) x 4(H) mm or 32(φ) x 14(H) mm	33(W) x 33(D) x 4(H) mm or 36(φ) x 14(H) mm
Supported SEM specimen mounts	10 mmφ x 2 pcs and 15 mmφ x 2 pcs and 32 mmφ x 1 pc.	10 mmφ x 7 pcs or 15 mmφ x 4 pcs or 36 mmφ x 1 pc.
Coating thickness	Several nanometers to hundreds of nanometers	
Settable min. thickness	Ultra-thin film mode: 0.1 nm / Normal mode: 1 nm	
Coating type	Osmium conductive film, Plasma-polymerized film (naphthalene)	Osmium film
OsO ₄ reservoir features	Detachable (capable to store in freezer) / Gas port integrated safety locking pin / Observation window / Built-in ampoule cutter	
Naphthalene reservoir features	Heater for effective sublimation / Observation window	-
Gas introduce/exhaust method	Automated control with vacuum gauge, solenoids, and vacuum pump	
Power requirements	100VAC (single-phase) 50/60Hz 12A (including supply to a vacuum pump)	100VAC (single-phase) 50/60Hz 10A (including supply to a vacuum pump)
Dimensions	450(w) x 410(D) x 390(H) mm	450(w) x 390(D) x 340(H) mm
Weight	Approx. 30kg	Approx. 20kg

Vacuum rotary pump (supplied with Osmium absorption filter)

Power requirements	100VAC (single-phase) 50Hz(60Hz) 550W	100VAC (single-phase) 50Hz(60Hz) 200W
Full-load current	9.0A (8.4A)	5.6A (4.8A)
Actual pumping speed	200L(240L) /min	50L(60L) /min
Dimensions	170(W) x 520(L) x 560(H) mm	170(W) x 400(L) x 580(H) mm
Weight	Approx. 31kg	Approx. 18kg

Specifications – optional units

Conductive ultra-thin film coating system (low-current method)

Coating type	Osmium conductive ultra-thin film	
Coating thickness	0.5~3.0 nm	
Gas introduce method	Needle valve	

Conductive ultra-thin film coating system (mixed-gas method)

Coating type	Osmium conductive ultra-thin film	-
Gas introduce method	Mass-flow controller	-
Supported gas	N ₂ , Ar, He, O ₂	-
Dimensions	570(w) x 410(D) x 390(H) mm	-

*This option cannot be installed along with hydrophilizing treatment system

Hydrophilizing treatment system

Gas introduce method	Needle valve	
Supported gas	Air	
Dimensions	570(w) x 410(D) x 390(H) mm	570(w) x 390(D) x 340(H) mm

*This option cannot be installed along with Conductive ultra-thin film coating system (mixed-gas method)

Deep well electrode


Reaction chamber	Glass chamber, 160(φ) x 125(H) mm	Glass chamber, 120(φ) x 103(H) mm
Max. sample size	44(W) x 44(D) x 4(H) mm or 46(φ) x 44(H) mm	33(W) x 33(D) x 4(H) mm or 36(φ) x 44(H) mm
Supported SEM specimen mounts	10 mmφ x 16 pcs or 15 mmφ x 7 pcs or 46 mmφ x 1 pc.	10 mmφ x 7 pcs or 15 mmφ x 4 pcs or 36 mmφ x 1 pc.
Dimensions	450(w) x 410(D) x 410(H) mm	450(w) x 390(D) x 370(H) mm

Technical data

Comparison data: Osmium plasma coat vs. Sputter coat

Sample: Human Lymphocyte


Filgen's OPC



Direct Magnification: 10,000X


Coating Device: OPC
Coating Material: Osmium
Coating Thickness: 5 nm
Accelerating Voltage: 15.0 kV

Highly smooth coating surface



Direct Magnification: 50,000X


Ion beam sputter coater



Direct Magnification: 10,000X

Coating Device: Ion Beam Sputter
Coating Material: Pt-Pd
Coating Thickness: 8 nm
Accelerating Voltage: 15.0 kV


Grains of Pt-Pd are observed at higher magnification




Direct Magnification: 50,000X

Sample: Barium Titanate

Filgen's OPC




Direct Magnification: 10,000X




Direct Magnification: 100,000X

Coating Device: Osmium Plasma Coater
Coating Material: Osmium
Coating Thickness: 3 nm
Accelerating Voltage: 10.0 kV

Magnetron sputter coater



Direct Magnification: 10,000X




Direct Magnification: 100,000X


Coating Device: Magnetron Sputter
Coating Material: Pt-Pd
Coating Thickness: 3 nm
Accelerating Voltage: 10.0 kV

Horizontal bright and dark bands are observed due to charging

Ion beam sputter coater





Direct Magnification: 10,000X



Direct Magnification: 100,000X

Coating Device: Ion Beam Sputter
Coating Material: W
Coating Thickness: 1.5 nm
Accelerating Voltage: 10.0 kV

Horizontal bright and dark bands are observed due to charging





Hybrid coating with osmium film and plasma-polymerized film

By coating with plasma-polymerized film of naphthalene before osmium coating, a hybrid film which enables to protect samples, reduce peeling of resin block, as well as prevent surface charging is obtained. This hybrid film is suitable to prepare SEM samples of non-conductive hydrophobic material like fluoro-resin.

Comparison data: Osmium single coat vs. hybrid coat


Sample: fluoro-resin film, 0.3 mm thickness



Evaluation for repeatability

Sample: epoxy block

- Experimental protocol
1. Coat osmium and naphthalene alternately on an epoxy block,
Naphthalene film: 10 nm, 6 layers
Osmium film: 1.0 nm, 5 layers
 2. Re-embed the sample by epoxy resin,
 3. Prepare cross sections with ultra microtome, then observe with TEM.



For more details, please check our website: <https://filgen.jp/Product/Sl/English/OPC/technical.html>



Filgen, Inc.
Scientific Instruments Dept.

1-1409 Jonoyama, Midori-ku, Nagoya 459-8011
JAPAN
TEL: +81 52-624-4388
Email: si-support@filgen.jp
<https://filgen.jp/>

(OCT 2019)