

Datasheet 2021/06/22

Product Name: *Anti-Heteroxylan [LM11] Antibody*

Catalog Number: *ELD017*

Size: *5mL*

Specifications:

Product Type: Antibody

Antigen: Heteroxylan

Isotype: IgM

Clonality: Monoclonal

Clone Name: LM11

Reactivity: Land plants

Immunogen: Neoglycoprotein (xylopentaose-BSA)

Species Immunized: Rat

Epitope: Beta-1,4-xylan

Buffer: Cell Culture Supernatant, 0.05% Sodium Azide

Tested Applications: IF (1:10), ELISA (1:10)

Storage: +4C (short term), -20C (long term)

Shipped: Ambient temperature

Provider: *From the laboratory of [Paul Knox, PhD](#), University of Leeds.*

References:

1. McCartney L, Marcus SE, Knox JP. Monoclonal antibodies to plant cell wall xylans and arabinoxylans. *J Histochem Cytochem.* 2005 Apr;53(4):543-6.

If you publish research with this product, please [let us know](#) so we can cite your paper.

FOR RESEARCH USE ONLY. NOT INTENDED FOR DIAGNOSTIC or THERAPEUTIC USE.