

Datasheet 2021/06/22

Product Name: *Anti-Arabinogalactan-protein (AGP) [LM2] Antibody*

Catalog Number: *ELD021*

Size: *5mL (supernatant)*

Specifications:

Product Type: Antibody

Antigen: Arabinogalactan-protein (AGP)

Isotype: IgM

Clonality: Monoclonal

Clone Name: LM2

Reactivity: Land plants

Immunogen: Polysaccharide

Species Immunized: Rat

Epitope: Glucuronosyl residue of AGP

Buffer: Cell Culture Supernatant, 0.05% Sodium Azide

Tested Applications: IF (1:10), ELISA (1:10)

Storage: +4C (short term), -20C (long term)

Shipped: Ambient temperature

Provider: *From the laboratory of [Paul Knox, PhD](#), University of Leeds.*

References:

1. *Smallwood M, Yates EA, Willats WGT, Martin H, Knox JP (1996) Immunochemical comparison of membrane-associated and secreted arabinogalactan-proteins in rice and carrot. Planta 198, 452-459*

If you publish research with this product, please [let us know](#) so we can cite your paper.

FOR RESEARCH USE ONLY. NOT INTENDED FOR DIAGNOSTIC or THERAPEUTIC USE.