

Datasheet 2021/06/22

Product Name: *Anti-Grass Heteroxylan [LM27] Antibody*

Catalog Number: *ELD034*

Size: *5mL (supernatant)*

Specifications:

Product Type: Antibody

Antigen: Heteroxylan

Isotype: IgM

Clonality: Monoclonal

Clone Name: LM27

Reactivity: Grass

Immunogen: Polysaccharide

Species Immunized: Rat

Buffer: Cell Culture Supernatant, 0.05% Sodium Azide

Tested Applications: IF 1:10

Storage: -20C

Shipped: Room Temperature

Provider: *From the laboratory of [Paul Knox, PhD](#), University of Leeds.*

References:

1. Cornuault V, Buffetto F, Rydahl MG, et al. Monoclonal antibodies indicate low-abundance links between heteroxylan and other glycans of plant cell walls. *Planta*. 2015;242(6):1321-1334.

If you publish research with this product, please [let us know](#) so we can cite your paper.

FOR RESEARCH USE ONLY. NOT INTENDED FOR DIAGNOSTIC or THERAPEUTIC USE.